



CENTRON
*AccuMix*TM **CM760**
DATA COLLECTION MIXER

Centron

Smart Solutions for Smarter Bloodbanking



Redefining the standards

AccuMix™ CM760

DATA COLLECTION MIXER

Features

- Accurate weighing and mixing of collected blood
- 5-inch LCD screen with user-friendly interface
- Visual/audible alarms with voice guide
- Multi-functional remote control unit
- Dual barcode scanner system
- Wireless data transfer (CDM software)
- Emergency onboard battery
- Large capacity battery pack for portable use
- *SureSeal™* tube sealer integration ready
- Secured stackable storage - Carrying Case Pro
- Versatile configuration



* shown in full configuration (CM760 Extended)

One for All, All for One

Multi-functional Remote Control Unit

- 10 LED lamps indicate % progress of collection and alarms. The lamps flash in case of alarming.
- Main operation buttons (CLAMP/START/END/ENTER) on the unit enhance user convenience.
- Secondary barcode scanner is equipped at the front.





SureSeal™ Tube Sealer Integration

- SureSeal™ tube sealer module can be installed in the main unit. By plugging in a sealing hand unit, CM760 is ready to offer the highest level of versatility.



High User Customizability

- CM760 is designed to satisfy various operating conditions in professional demands. Users can adjust most operational setting including target collection volumes, alarm levels and on/off, tray moving speed and even specific gravity of blood.
- Centron R&D Lab offers customization support to meet your specific operation requirements.



Intelligent Voice Guide

- Voice guide system speaks real-time information including selected target volume, collection start, collection end and alarm conditions.
- 10 selectable languages:
 - English, Spanish, French, German, Portuguese, Italian, Russian, Chinese, Korean and Japanese



Dual Barcode Scanner System

- Barcode scanners read most of standard barcode labels including blood bank-specified ISBT128 code.
- Centron R&D Lab is ready for customized solutions.



Dual Battery System

- Onboard battery (rechargeable) is for emergency backup. In case of unexpected AC power failure, it guarantees saving all donor data and supplies power for 30 more minutes to complete ongoing collection cycle. No need to use BN128 battery pack all the time in case of AC powered indoor uses.
- For portable or outdoor use, the large capacity BN128 battery pack (rechargeable & replaceable) is available. The pack supports up to 8 hours of operation without AC mains.
- BN128 battery pack is charged through the main unit or BC164 fast battery charger.



■ 5-inch LCD Screen

- Clear backlit LCD display helps users following instructions and monitoring data. Display contains all necessary information for user convenience.




■ User-friendly Interface

- Collection processes are represented with graphic symbols.
- Users can change settings depending upon requirements with passwords.




During Collection

- Collected blood volume
- Target blood volume
- Total weight
- Collection elapsed time
- Blood flow rate and alarm
- Collection progress bar



Collected volume curve

- presents collected blood volume vs elapsed time
- Curve slope gives a quick view of collection flow rate




Data summary and review


- shows donor data summary at the end of collection cycle
- scrolls all data saved in the internal device memory

■ Smart Data Management System

- All the collection data including up to 1,000 donor information are saved in the internal device memory of CM760 main unit. Each data consists of 30 information fields.
- Data can be transferred by a USB memory stick or via wireless data link (CDM).
- One CDM network can group 18 units. A supervisor can monitor and control real-time operation of up to 18 units. Donor data, collection data, alarm status, device settings are real-time managed by a computer.
- Data format is available in .txt (text file) and .csv (general spreadsheet file). Users can manage data with a spreadsheet program. Microsoft Excel® is recommended.
- Device settings in the network can be managed as all or separately.



| No. | Donor ID | Device ID | Status | Preset | Collected | Progress | Time |
|-----|------------|------------|------------|--------|-----------|----------|-------|
| 1 | [REDACTED] | [REDACTED] | Ready | 0 ml | 0 ml | 0 % | 00:00 |
| 2 | [REDACTED] | G760ICM42 | Collecting | 450 ml | 50 ml | 11 % | 01:23 |
| 3 | [REDACTED] | G760ICM43 | Low Flow | 450 ml | 220 ml | 51 % | 06:50 |
| 4 | [REDACTED] | G760ICM44 | Collecting | 450 ml | 234 ml | 50 % | 03:42 |
| 5 | [REDACTED] | G760ICM45 | Ready | | | | |
| 6 | | | | | | | |
| 7 | [REDACTED] | G760ICM47 | End | 450 ml | 450 ml | 100 % | 08:46 |
| 8 | [REDACTED] | G760ICM48 | Ready | | | | |
| 9 | [REDACTED] | G760ICM49 | High Flow | 450 ml | 300 ml | 66 % | 07:05 |
| 10 | | | | | | | |
| 11 | [REDACTED] | G760ICM11 | T1 Pass | 450 ml | 420 ml | 93 % | 10:02 |
| 12 | [REDACTED] | | | | | | |
| 13 | | | | | | | |
| 14 | | | | | | | |
| 15 | | | | | | | |
| 16 | | | | | | | |
| 17 | | | | | | | |
| 18 | | | | | | | |



CENTRON
AccuMix™ CM760
 DATA COLLECTION MIXER

Optional Chart

| | Selective Item | Including | Note |
|--------------------|--|--|---|
| Unit Configuration | • CM760 Basic | • CM760 main unit • AC power adapter • User's manual | • Essential configuration for BASIC MODE operation ONLY • BN128 battery pack is not included |
| | • CM760 Standard | • CM760 Basic • [Option D] included • [Option C] included | • Standard configuration for DATA MODE operation • BN128 battery pack is not included |
| | • CM760 Extended | • CM760 Standard • [Option S] included • BN128 battery pack | • Complete configuration with tube sealing feature |
| Optional Package | • [Option D] DATA Management Package | • Barcode scanner • Multi-functional remote control unit • USB memory stick | • Supports DATA MODE operation • Either [Option C] or Pole Mounting Kit may be required to hang a barcode scanner, multi-functional remote control unit and sealing hand unit |
| | • [Option C] Carrying Case Package | • Carrying Case Pro (including device mounting platform, pole, hanger & thumb screw, plug extension adapter) | • Stackable and secured storage • Easy transportation of all parts • Main unit charging available with the case closed • It can be used as collection stand for main unit • Pole and hanger are mounted on device mounting platform |
| | • [Option S] Tube Sealer Package | • Built-in <i>SureSeal™</i> tube sealer module • Sealing hand unit & 1.8m(6ft) cable | • Built-in <i>SureSeal™</i> tube sealer module (SEM 01) • BN128 battery pack must be installed |
| Supporting Device | • Collection Data Manager (CDM) | • CDM module • Antenna & USB cable • Installation CD | • Wireless Modem with CDM software package • Real-time collection site monitoring for up to 18 units • Wireless data transferring to a computer • Simultaneous device setting for a group of 18 units |
| | • BC164 Fast Battery Charger | | • AC powered fast multi-battery charger • Charging of 6 BN128 battery packs |

* Optional parts including BN128 battery pack may be sold separately. For more information about package customization, please visit our website.

Specification

| | |
|-----------|--|
| Product | Data Collection Mixer |
| Model | CM760 |
| Power | 18V--- 1.67A + - |
| Adapter | 100-240V~ 50/60Hz 1.5A |
| Dimension | Main Unit: 200(W)mm/400(L)mm/176(H)mm Carrying Case: 280(W)mm/450(L)mm/380(H)mm |

Information



DISTRIBUTOR

Centron Centron Corporation

Address: 44 Beolmal-ro, Dongan-gu, Anyang, Gyeonggi-do,
Postal Code 14059 / Republic of Korea
Phone: +82-31-422-7890 / **Fax:** +82-31-422-7866
E-mail: info@centrons.com
Website: www.centrons.com



* The specification and design are subject to change without prior notice for product upgrades.