


Don't just seal, SURESEAL



CENTRON

SureSeal™ SE730

TUBE SEALER

Pack-n-Go

SureSeal™ SE730 is a compact portable sealer

■ Features

New Li-Ion New Li-Ion Battery - Fast rechargeable and replaceable battery managed by automatic sealing output control software


- The world first customizable adjustable sealing power modes
- Small and lightweight Li-Ion battery can make up to 800 seals
- Seals at *SureSeal™* hand unit
- Seal wide range of blood bag tubing - optimal for tubing diameter 2-6 mm (1/4")
- No warm-up required / Fast sealing within one second
- Detachable head protector for easy cleaning
- All inclusive durable carrying bag - ideal for mobile use
- Built-in battery charger

Centron

Smart Solutions for Smarter Bloodbanking

Control panel

- Indication lamps on the front panel indicate available battery charge and charging status. For stable and continuous operation, it is recommended that the battery pack should be recharged when two or more indication lamps are lighting. Recharging to full capacity takes 3 hours.



Battery charge level lamps

- The five green LED lamps indicate battery charge level stored in the battery pack inside. Each lamp denotes 20 % of full charge capacity of the battery.

Charging status lamp

- The yellow LED lamp turns on when the charger is plugged into the charging jack. The lamp flashes when charge level reaches to full.

Sealing lamp

- The blue LED lamp indicates RF power is being applied to the hand unit. The lamp turns on by pulling the trigger, and turns off automatically when a seal is made, or after 5 seconds.

Sealing hand unit & Carrying bag

- Additional plug-in type *SureSeal™* sealing hand unit
- Ergonomic design with standard 1.8 m (6 ft) cable
- Satisfying various operating conditions
- Carrying bag houses all parts



Specification

Product	TUBE SEALER
Model	SE730
Power	18V --- 1A + -
Battery	Li-Ion 14.4V/2.5AH
Adapter	100-240V~ 50/60Hz 0.6A
Dimension	Main unit: 152(W)mm/194(L)mm/39(H)mm Hand unit: 28(W)mm/204(L)mm/44(H)mm

Information



DISTRIBUTOR

Centron Centron Corporation

Address: 44 Beolmal-ro, Dongan-gu, Anyang, Gyeonggi-do,
Postal Code 14059 / Republic of Korea
Phone: +82-31-422-7890 / **Fax:** +82-31-422-7866
E-mail: info@centrons.com
Website: www.centrons.com

